HELLO KITTY

ETB-134956.1

Elektrische Zahnbürste für Kinder (D)
Elektrische tandenborstel voor kinderen (NL)
Cepillo de dientes eléctrico para niños (E)
Escova de dentes elétrica para criança (P)
Spazzolino elettrico per bambini (I)
Detská elektrická zubná kefka (SK)
Otroška električna zobna ščetka (SLO)
Elektryczna szczoteczka do zębów dla dzieci (PL)
Periuță de dinți electrică pentru copii (RO)
Dječja električna četkica za zube (HR)
Elektrický zubní kartáček pro děti (CZ)
Електрическа четка за зъби за деца (BG)
Gyermek elektromos fogkefe (H)
Kids electric toothbrush (EN)

Inhalt – Inhoud – Contenido – Índice – Contenuto – Obsah – Vsebina –
 Treść – Cuprins – Sadržaj – Obsah – Съдържание – Tartalom – Content

Bedienungsanleitung – German (D)	- 2 -
Gebruiksaanwijzing – Netherlands (NL)	- 9 -
Manual de Instrucciones – Spanish (E).....	- 16 -
Manual de instruções – Portuguese (P)	- 23 -
Manuale d'istruzioni – Italian (I).....	- 30 -
Návod na obsluhu – Slovak (SK)	- 36 -
Navodila za uporabo – Slovenian (SLO).....	- 42 -
Instrukcja obsługi – Polish (PL)	- 48 -
Manual de utilizare – Romanian (RO)	- 56 -
Uputa za korištenje – Croatian (HR)	- 62 -
Návod k použití – Czech (CZ)	- 68 -
Ръководство с инструкции – Bulgary (BG)	- 74 -
Kezelési útmutató – Hungary (H)	- 81 -
Instruction manual – English (EN)	- 87 -

SICHERHEITSINFORMATIONEN

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SAFETY INSTRUCTIONS

Before use make sure to read all of the below instructions in order to avoid injury or damage, and to get the best results from the appliance. Make sure to keep this manual in a safe place. If you give or transfer this appliance to someone else make sure to also include this manual.


In case of damage caused by user failing to follow the instructions in this manual the warranty will be void. The manufacturer/importer accepts no liability for damages caused by failure to follow the manual, a negligent use or use not in accordance with the requirements of this manual.

1. Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision.
2. Children shall not play with the appliance.
3. Check and make sure to use an adapter with a rated output voltage of 5V, a rated output current of 1A, and a waterproof level of IPX4 to ensure the electrical safety. The adapter cannot be used for charging this product if it does not comply with the above requirements.
4. If the USB supply cord with charging base is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
5. This appliance contains batteries that are non-replaceable.
6. If the appliance drops into water while it is being charged, please firstly disconnect the power source and then remove the appliance from the water.
7. Please charge the appliance in a dry environment far away from water. Do not charge the appliance in a humid environment such as bathroom and others.
8. It is forbidden to use and rinse the toothbrush while it is being charged.
9. Keep your hands dry when handling the power cable for charging.
10. Do not place the appliance and its battery in a microwave oven or on an induction cooker.
11. Do not disassemble or modify the toothbrush before disposing of it for recycling.
12. Do not immerse the toothbrush into hot water.
13. Do not press against tooth or gum with the toothbrush using excessive force. Do not brush the same position for long periods time to avoid any damage.
14. **WARNING:** Make sure that there is no appearance of damage before each use. Do not use the toothbrush if damaged.
15. **WARNING:** Please consult your doctor whether you can use the toothbrush If you have had dentistry surgical operation in past two months or are under being subject to dental treatment.
16. **WARNING:** Oral hemorrhage may occur when using this product at the first time. If the oral hemorrhage lasts for more than a week, please consult your doctor.
17. **WARNING:** Please consult your doctor before using this product if you are using a heart pacemaker or other electronic devices on or in your body.
18. **WARNING:** Please consult your doctor before using this product if you have medical problems.
19. **HAZARD:** Do not heat the toothbrush, throw it into the fire, store, charge or use at a high-temperature environment to avoid overheating, fire or explosion as the handle of this product is equipped with a built-in rechargeable battery.
20. The toothbrush is for personal use, do not share it with other persons.
21. This product is designed for cleaning teeth, gums, and tongue only.
22. The temperature for charging, using, and storing the product shall be between 5°C and 35°C.
23. Please recharge this product in time when the battery is low.
24. Do not use the toothbrush when you are taking a bath or shower.
25. This appliance must only be supplied at safety extra low voltage corresponding to the marking on the appliance.

26. During use and charging, the device may become warm, which is normal.
27. Do not place near heat sources such as radiators, cookers, etc.
28. Remove any oral jewellery prior to use of this appliance.
29. Do not attempt to repair, disassemble or alter the appliance yourself.
30. Do not use the appliance if you have sores, blisters or cuts in your mouth.
31. Regarding the information pertaining to the use of the appliance and user maintenance, please refer to the below paragraph of the manual.

All images shown in this manual are for reference purposes only and may vary from the actual product.

PARTS DESCRIPTION

1. Toothbrush head
2. On/off button
3. Mode button
4. Whitening mode indicator light
5. Clean mode indicator light
6. Sensitive mode indicator light
7. Battery indicator light
8. Handle
9. Charging base
10. USB cable

BEFORE USE

1. Charging your toothbrush with an adapter

- 1) Turn off the toothbrush before charging. Put it onto the charging base.
- 2) Place the charging base on a flat surface. Connect the USB cable into an adapter. Plug the adapter into an appropriate power source.

Note: Check and make sure to use an adapter with a rated output voltage of 5V, a rated output current of 1A, and a waterproof level of IPX4 to ensure the electrical safety. The adapter cannot be used for charging this product if it does not comply with the above requirements.

- 3) During charging, the battery indicator light will illuminate red. Once the toothbrush is fully charged, the indicator light turns green.

Note: The indicator light will flash red when the battery is low, which means that the toothbrush needs to be recharged.

Note: It may take approx. 12-15 hours to fully charge the toothbrush.

Notes:

- When the toothbrush is charged for the first time or if the toothbrush has not been used for a long time, the indicator light may not illuminate for the first 5 minutes during charging. This is caused by the low battery of the toothbrush, which is normal.
- The toothbrush must be charged completely under the following conditions:
 - The toothbrush is used for the first time.
 - The toothbrush has not been used for more than 3 months.
 - The toothbrush vibration becomes weak or low-battery indicator light is shown.
- The toothbrush may become hot while being charged. This is normal and will not affect the performance of the toothbrush.

Low-battery reminding function:

When the battery becomes low during operation, the battery indicator light will flash continuously until the unit powers off.

2. Installing the toothbrush head

Push the toothbrush head downwards onto the handle. Make sure they are properly connected. To remove the toothbrush head, pull it upwards.

USE

1. Wet the brush head and apply some toothpaste.
2. Place the toothbrush at a 45-degree angle against the gumline of your teeth.
3. Press the on/off button to turn on the toothbrush. Press the mode button to choose your desired cleaning mode. The toothbrush will operate as follows: Whitening mode, use with whitening toothpaste – Normal clean mode for daily brushing – Sensitive mode for sensitive teeth.
4. Gently move the toothbrush head back and forth between the teeth slowly. The toothbrush will pause for a short while approx. every 30 seconds to remind you to brush other sections of your mouth. The toothbrush will keep running at the current cleaning mode and shut off automatically in 2 minutes. The current cleaning mode will be saved and the toothbrush will work according to the saved cleaning mode when turned on again.
5. Once finish brushing your teeth, turn off the toothbrush. Rinse the brush head and bristles under running water, then leave them to dry. Wipe the handle with a towel.

Note: Do not scrub hard. Apply light pressure to maximize the toothbrush's effectiveness.

CLEANING AND MAINTENANCE

CAUTION: Before cleaning, disconnect the charging base from the power and remove the toothbrush.

Rinse the toothbrush head with clean water at a temperature lower than 45°C and then dry the surface of the toothbrush handle with a towel.

CAUTION: Do not immerse the toothbrush handle into water. Never use poisonous cleaning agents.

CAUTION: Do not put any part of the toothbrush into the dishwasher for cleaning.

CAUTION: Do not use chemical, steel, wooden or abrasive cleaners to clean the toothbrush.

CAUTION: Please do not use a deformed toothbrush head or a toothbrush head with bent bristles.

CAUTION: Replace the toothbrush head every 3 months or sooner if signs of wear appear.

Cleaning the toothbrush head

1. Rinse the brush head and bristles under running water after each use.
2. Detach the brush head and clean its inner surface with warm water at least once a week.
3. Dry thoroughly.

Cleaning the handle

1. Detach the toothbrush head and rinse the metal shaft with warm water. Make sure all the residual toothpaste is removed.
2. Wipe the handle surface with a damp cloth and then dry thoroughly.

Storage

If the toothbrush will not be used for a long time, disconnect the power source and clean the toothbrush accordingly. Store the toothbrush in a cool and dry place away from the direct sunlight. Charge the toothbrush every three months to prevent the battery from discharging completely. A fully-discharged battery may not be able to be recharged.

TROUBLESHOOTING

Problem	Cause	Solution
Reduced cleaning intensity	Excessive pressure on the teeth	Gently place the brush head on the teeth.
The toothbrush works for only a short time even after charging.	1. Charging time is insufficient. 2. The product has not been used for more than 6 months. 3. The battery has reached its lifespan.	1. Recharge the appliance to full power. 2. Contact customer service.
The toothbrush stops working.	The battery has run out.	Please recharge the battery.

TECHNICAL DATA

Main unit:

Rated operating power: Max 2W

Rated operating voltage: DC 3.7V

Rated input: 5V  0.1A

Charging base:

Rated charging power: 5W

Rated input: 5V  1A

Battery voltage: 3.7Vdc

Battery capacity: 600mAh

Use only with supplied charging base

IPX7

Working time: up to 18 days (2 times/day; 2 mins/time)

Replacement brushes:

Part number of Hello Kitty brushes:

Part-134980

Compatible with Oral-B Pro® brush heads.

Oral-B Pro® is a registered trademark of its respective owner.

This product is not manufactured or endorsed by Oral-B.

Your EU legal guarantee rights

You have a legal guarantee for at least 3 years from delivery. If the product is faulty or not as described, the seller you bought from must provide a free-of-charge remedy: repair or replacement, or—if not possible within a reasonable time and without significant inconvenience—a price reduction or refund. This does not affect your statutory rights under EU and national laws.

Note: If a defect appears within 1 year of delivery, it is presumed to have existed at delivery unless the seller proves otherwise (in some countries this period is 2 years).

Exclusions/limitations

This guarantee does not cover consumables or parts subject to normal wear and tear unless the defect is due to lack of conformity. Physical damage caused by misuse or by not following the instruction manual is not covered. Nothing in this guarantee excludes or limits liability where such exclusion or limitation would be unlawful.

How to claim

Contact your seller/retailer with proof of purchase (receipt, invoice, order confirmation or other reasonable evidence).

Remedies under the legal guarantee are free of charge, including shipping, labour, and materials.

Electrical and electronic devices/ information for private households

The Directive 2012/19/EU on waste electrical and electronic equipment contains a large number of requirements for handling electrical and electronic equipment. The most important are compiled here.

1. Separate collection of waste electrical and electronic equipment

Electrical and electronic equipment that has become waste is referred to as waste electrical and electronic equipment (WEEE). Owners of WEEE must collect them separately from unsorted municipal waste. WEEE do not belong in household waste, but in special collection and return systems.

2. Batteries, accumulators and lamps

Owners of WEEE must hand over them in complete condition. However, owners of WEEE have to separate old batteries and accumulators which are not enclosed by the WEEE, as well as lamps which can be removed, from the WEEE without destroying it before handing it over to a collection point.

3. Ways to return the WEEE

Owners of WEEE from private households can hand them in free of charge at the designated collection points of the public waste disposal authorities or at the take-back points set up by manufacturers or distributors in terms of the Directive 2012/19/EU.

4. Privacy Protection Notice

WEEE often contain sensitive personal data. This applies especially to information and telecommunications technology devices such as computers and smartphones. In your own interest, please note that each end user is responsible for deleting the data on the WEEE before disposing it.

5. Potential environmental impact

WEEE contain substances that may have a negative impact on the environment and human health if their collection, hand-over, re-use or material recovery is not done in compliance with respective legislation.

6. Your role in WEEE waste treatment

By fulfilling these obligations and especially by fulfilling the obligation to collect WEEE separately, not to dispose it together with unsorted municipal waste and to hand it over to at collection points, you as the end user contribute to the re-use and material recovery of WEEE.

7. Meaning of the crossed-out dustbin symbol

The crossed-out dustbin, which is regularly displayed on electrical and electronic equipment, indicates that the respective device must be collected separately from unsorted municipal waste at the end of its life cycle.

BATTERY MUST BE RECYCLED OR DISPOSED OF PROPERLY. DO NOT OPEN. DO NOT DISPOSE OR IN FIRE OR SHORT CIRCUIT.

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KiK Textilien und Non-Food GmbH
Emerio B.V.