bebeconfort

SYSTECH HEALTHCARE CO., LTD.
Room 1901, Huamin Empire Plaza,
Economic Development Zone, 311100,
Hangzhou, Zhejiang, China

EC REP Shanghai International Holding Corp. (Europe) GmbH
Eiffestraße 80, 20537 Hamburg, Germany

Imported and Imported by STIL - Cité Pariwest, 7 7000 Vaux-le Penil - France

CE 0123 REF 3106203000
3106203100
3106203101

Keep dry. A conserver dans un endroit sec. Bleib trocken. Droog houden.
Conserve in luogo asciutto. Mantener en lugar seco. Manter seco. Udržujte v suchu. Uchowywać w suchu. Ne déposez pas de pression. Не давите.

LCD On/Off button

Probe

Battery cap

V2 - 06/09/2024

EN - ULTRA-FAST FLEXIBLE THERMOMETER

PLEASE READ CAREFULLY BEFORE USING :
- This digital thermometer provides a quick and highly accurate reading of an individual's body temperature. The digital thermometer is intended to measure the human body's temperature through oral, rectal or axillary (under the arm) use. It is suitable for use on people of all ages, including children under 3 years old. Parents or caregivers are responsible for using this device and for setting the correct mode for the person whose temperature is being taken.
- ISO 80601-2-56:2017 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement. ASTM E1112-00(2013) Standard Specification for Electronic Thermometer for Intermittent Determination and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of IEC 60601-1-11:2015+AMD1:2020+AMD2:2020 (Edition2.0)Safety standards. And the manufacturer is ISO 13485 certified.

Contents : 1 thermometer, 1 owner's manual and 1 storage case.

- The flexible probe thermometer is designed to be used exclusively to measure body temperature in the human body. It has a high sensitivity and a high level of accuracy. Its flexible probe remains the most suitable for the use of small children.
- The thermometer is designed to operate within the following ranges: Temperature of ages, Humidity and Atmospheric pressure.

Please carefully read this guide before using the thermometer and keep it in a safe place for reference.

SYMBOL EXPLANATION

Caution!		LOT	Batch Code
Consult Accompanying Documents			Manufacturer
MD	Medical Device	-10°C 70°C	Storage and Transportation Temperature Limit: -20°C-55°C (-4°F ~ 131°F)
CE 0123	European Authorized Organization		Storage and Transportation Humidity Limit: 15%~95% RH
	Disposal of this product and used batteries should be carried out in accordance with the relevant regulations for the disposal of electronic products.		
CE	The product is in compliance with the requirements of the regulation for medical devices (EU) 2017/745. 0123 is the identification number of the notified body.		
IP27	IP27 - Number 2 : Protected against solid foreign objects of 12.5mm Ø and greater IP - Number 7 : Protected against the effects of temporary immersion in water		

A. TEMPERATURE MEASURING METHODS

1. Rectal temperature

Lubricate the tip with Vaseline or water to facilitate insertion. Gently insert the thermometer about 2.5cm (1 inch) into the rectum. Hold it in place for the duration of the reading.

2. Axillar temperature

Wipe armpit with dry towel. Move the tip under the armpit into the hollow and press arm down firmly. Hold it in place for the duration of the reading. Keep arm down during reading.

3. Oral temperature (Avoid drinking hot or cold beverages)

Place its head under tongue against frenum. Recline mouth and breath very slowly through nose. Hold it in place for the duration of the reading. Keep mouth closed.

B. CONTAMINATION

Do not use in case of major lesions affecting visual measuring areas, on a very nervous patient, in case of risk of ingestion or inhalation.

C. CLEANING AND DISINFECTION

Before and after each usage, clean the tip with warm soapy water or a soft cleaning product. Rinse with clean water and dry with a soft cloth. Do not use any abrasive or harsh cleaning solutions. The thermometer must always be in good condition (no cracks or deformation) and into a disinfectant liquid.

D. TEMPERATURE READING

Press the On/Off button (1 beep) immediately. The thermometer will then switch to °C mode. Then the last recorded temperature is shown after it is first used and then [88.8] the flow disappeared and the °C indicator is shown. Wait for the thermometer to reach the temperature of where it's the maximum output power rating of the transmitter in watts (w) according to the transmitter classification. Recommended separation distance between portable and mobile RF communications equipment and the medical device. Ensure that children do not use the instrument unsupervised. Some parts are small enough to be swallowed. Use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 1.5m from such devices when using this device. • Remove battery from the device when not in operation for a long time.

E. GOOD OPERATION CONDITIONS

- Ensure that the ambient temperature is stable before using the thermometer.

1. Display field test:

- Quick press and release [ON/OFF] to power: display field test is performed.

2. Measuring system test:

- Press and release [ON/OFF] button immediately then [88.8] with the °C flashing after the thermometer switches to measurement system. However, if [88.8] screen displays [Lo] the thermometer cannot be used and should be replaced.

- H_2 Temperature taken is higher than 42.9°C(109.9°F). Turn off, wait one minute and take a new temperature via close contact and sufficient rest.

Error message	Problem	Solution
Lo	Temperature taken is lower than 32.0°C (89.6°F)	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
Hi	Temperature taken is higher than 42.9°C (109.9°F)	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
Err	The system is not working properly	Remove the battery, wait for 1 minute and replenish it. If the message occurs, contact the reliable for service.
	Empty battery. Battery icon is flashing. Battery level can't be measured.	Replace battery.

H. BATTERY REPLACEMENT

- Carefully remove the thermometer end cover.
- Replace with a new battery. Keep the battery out of children's object, do not leave within the reach of infants.
- Replace with a new battery, as indicated by a low battery type SR 44, LR 44, or LR 1/2 357, 1.5v positive side (+) facing up.
- Replace the cover.
- Dispose of the old battery or the used thermometer in an appropriate manner.

I. TECHNICAL SPECIFICATIONS

- Measuring range: 32.0°C-42.9°C / 89.6°F-109.4°F ambient operating range
- Accuracy: ± 0.1°C for other measuring and ambient operating range
- Scale: 0.1°C or °F
- Operating mode
- Power supply: ≤ 1.50 V button battery
- Battery life: approximately 2000 measurements (using new batteries) ⚠
- Signaling time: 1 beep during turn on and off. 10 beeps at end of measurement.
- Signaling tone: 3 beep during turn on and off. 10 beeps at end of measurement.
- Memory recall: 1 value
- Ambient operating range: Temperature: 5°C-40°C; Relative humidity: 15%~95%RH; Atmospheric Pressure: 700hPa ~ 1060hPa
- Storage and transport conditions: Temperature: -20°C-55°C; Relative humidity: 15%~95%RH; Atmospheric Pressure : 700hPa ~ 1060hPa
- Weight: 15g
- Recommended service life: 3 years

J. WARRANTY

Your thermometer is guaranteed for a two-year period from purchase date. However, this warranty does not apply to damages, or any costs resulting from non-compliance with the guide, incorrect handling, misuse, shock or inappropriate occupation of the battery. It does not apply either to the battery or the package.

K. ELECTROMAGNETIC COMPATIBILITY

This device complies with the requirements of IEC 60601-1-2:2014+AMD1:2020 Edition 4.1 2020-12. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical device that requires special precautions regarding electromagnetic fields. It must be used only for the purpose for which it was intended. Portable and mobile RF communications equipment can affect the device. Use of the unit in conjunction with non- approved accessories can result in increased emissions or decreased immunity. Do not use adjacent to other devices. The device is not directly adjacent to or stacked on other devices.

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

L. CALIBRATION

The thermometer is calibrated during the manufacturing process. If used as per the instructions, it should not require calibration. However, if you suspect inaccurate readings, we suggest you check it every or whenever there is doubt about its clinical accuracy. To do so, please send the complete unit to the Customer Service Center. Please note that the warranty does not cover calibration which does not exist for regular requirements. You must always comply with the requirements related to the correct use of the instrument, functionality, and accuracy of the device as mandated by relevant laws and directives.

Table 1

Guidance and manufacturer's declaration - electromagnetic emission

Emissions test	Compliance	Electromagnetic environment - guidance
Emissions test CISPR 11	Group 1	The thermometer uses only RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Emissions test CISPR 11	Class B	The thermometer is suitable for use in all establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Table 2

Guidance and manufacturer's declaration - electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV ±4 kV ±8 kV ±15 kV air	±8 kV contact ±2 kV ±4 kV ±8 kV ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50Hz or 60Hz	30 A/m, 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3

Guidance and manufacturer's declaration - electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz to 2.70Hz	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz to 2.70Hz	Electromagnetic environment - guidance Portable and mobile RF communications equipment should be used no closer to any part of the thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. **Recommended separation distance** $d = \frac{3.5}{V_1}\sqrt{P}$ $d = \frac{3.5}{E_1}\sqrt{P}$ 80 MHz to 800 MHz $d = \frac{7}{E_1}\sqrt{P}$ 800 MHz to 2.7 MHz where P is the maximum output power rating of the transmitter in (w) according to the transmitter classification. Recommended separation distance between portable and mobile RF communications equipment and the medical device. Ensure that children do not use the instrument unsupervised. Some parts are small enough to be swallowed. Use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 1.5m from such devices when using this device. • Remove battery from the device when not in operation for a long time.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of communications equipment.

Rating of max. output of transmitter	Separation distance according to frequency of transmitter - m	
	80 MHz to 800 MHz $d = \frac{3.5}{V_1}\sqrt{P}$	800 MHz to 2.7 GHz $d = \frac{7}{E_1}\sqrt{P}$
0.01	0.04	0.07
0.1	0.12	0.23
1	0.35	0.7
10	1.11	2.22
100	3.5	7.0

For transmitters rated at a maximum output power not listed above the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer. NOTE 1 : At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

H_1 Temperature taken is higher than 42.9°C(109.9°F). Turn off, wait one minute and take a new temperature via close contact and sufficient rest.

H_2 Temperature taken is higher than 42.9°C(109.9°F). Turn off, wait one minute and take a new temperature via close contact and sufficient rest.

Err Error. The thermometer must be replaced. This is also true when the display field exists light flashes or does not turn on anymore.

Table 5

Recommended separation distances between RF wireless communications equipment

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of communications equipment.

Frequency	Max. Power	Distance	IEC 60601 Test Level	Compliance Level	Electro-magnetic environment - guidance
MHz	W	m	V/m	V/m	
385	1.8	0.3	27	27	RF wireless communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. **Recommended separation distance** $d = \frac{3.5}{E_1}\sqrt{P}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and E is the test level in V/m. Interference may occur in the vicinity of equipment marked with the following symbol:
450	2	0.3	28	28	
710					
745	0.2	0.3	9	9	
780					
810					
870	0.3	0.3	28	28	
930					
1720					
1845	0.3	0.3	28	28	
1970					
2450	2	0.3	28	28	
5240	0.2	0.3	9	9	
5500	0.2	0.3	9	9	
5785					

NOTE 1 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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